



## Traceable Calibration Services

Ensure your equipment operates successfully with fully-traceable Equipment Calibration Services.

Calibration is an essential part of equipment maintenance. STERIS offers calibration procedures to address your equipment's specific needs over its entire lifecycle.

Performed by specially trained, field-based technicians, all calibration tests are fully traceable to international standards and detailed in a complete report kept on file with STERIS for a minimum of three years.



**NIST**  
National Institute of  
Standards and Technology

Calibration and verification points are determined by STERIS, the manufacturer or your SOPs. We will make necessary adjustments to 'out of tolerance' conditions we may find and include the results in your final report.

**A full report is provided after every calibration. Your STERIS technician will review with you this report, which will contain:**



**A LIST**

of 'As Found' and 'Corrected' readings



**TOLERANCE PAGE AND SIGNATURE PAGE**

including a list of standards, procedures, temperatures, pressures, comments and signatures from both you and your STERIS technician



**CALIBRATION CERTIFICATES**

for each instrument used in the calibration

Your equipment will also be tagged to show cGMP compliance.



**Calibration is Offered for:**

- Sterilizers
- Washers
- Generators
- Water Stills
- Low temperature sterilization and biodecontamination (VHP) equipment

For more information on our Equipment Service options, please visit us at:  
[sterislifesciences.com/services/specialty-equipment-services-and-upgrades](https://sterislifesciences.com/services/specialty-equipment-services-and-upgrades)